

(12) INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(19) World Intellectual Property  
Organization  
International Bureau



(43) International Publication Date  
5 August 2004 (05.08.2004)

PCT

(10) International Publication Number  
**WO 2004/065929 A2**

(51) International Patent Classification<sup>7</sup>: **G01N**  
(21) International Application Number:  
PCT/EP2004/000313  
(22) International Filing Date: 16 January 2004 (16.01.2004)  
(25) Filing Language: English  
(26) Publication Language: English  
(30) Priority Data:  
03001057.3 17 January 2003 (17.01.2003) EP

(71) Applicant (for all designated States except US): EU-  
ROPÄISCHES LABORATORIUM FÜR MOLEKU-  
LAR BIOLOGIE (EMBL) [DE/DE]; Meyerhofstr. 1,  
69117 Heidelberg (DE).

(72) Inventor; and  
(75) Inventor/Applicant (for US only): FRANZ, Thomas  
[DE/]; Alte Eppelheimer Strasse 35 - 1, 69 115 Heidelberg  
(DE).

(74) Agents: WEICKMANN & WEICKMANN et al.; Post-  
fach 860 820, 81635 München (DE).

(81) Designated States (unless otherwise indicated, for every  
kind of national protection available): AE, AG, AL, AM,  
AT, AU, AZ, BA, BB, BG, BR, BW, BY, BZ, CA, CH, CN,  
CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, EG, ES, FI,  
GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE,  
KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD,  
MG, MK, MN, MW, MX, MZ, NA, NI, NO, NZ, OM, PG,  
PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, SY, TJ, TM,  
TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM,  
ZW.

(84) Designated States (unless otherwise indicated, for every  
kind of regional protection available): ARIPO (BW, GH,  
GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW),  
Eurasian (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), Euro-  
pean (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FI, FR,  
GB, GR, HU, IE, IT, LU, MC, NL, PT, RO, SE, SI, SK,  
TR), OAPI (BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW,  
ML, MR, NE, SN, TD, TG).

Published:  
— without international search report and to be republished  
upon receipt of that report

For two-letter codes and other abbreviations, refer to the "Guid-  
ance Notes on Codes and Abbreviations" appearing at the begin-  
ning of each regular issue of the PCT Gazette.

(54) Title: SILICONE/GRAPHITE SAMPLE HOLDER

(57) Abstract: The present invention relates to a sample holder for a mass spectrometer onto which a mixture of silicone and graphite is applied.



WO 2004/065929 A2

**Silicone/graphite sample holder****Description**

5

The present invention relates to a sample holder for a mass spectrometer onto which a mixture of silicone and graphite is applied.

10

15

20

25

The analysis of proteins by mass spectrometry has become a standard procedure in molecular biology in recent years. Peptide samples can be prepared by digestion of purified proteins directly or by an in-gel digestion of proteins previously separated by 1D or 2D-gel electrophoresis and mixed with a matrix for further analysis. A holder (target) is used onto which the samples are spotted. Presently, there are several robots available for spotting but for most applications manual spotting is necessary. Several types of holders exist. The most common one is a steel holder with ring shaped grooves. The samples are spotted into the ring, which is meant to prevent the samples from leaking out and cross-contaminating each other. Therefore, the volume that can be applied to that type of holder is fairly small and the crystals are very far apart. Additionally, the peptides bind only weakly to the steel surface and therefore the samples cannot be washed to get rid of surplus salt after spotting. This often results in poor data in the following mass spectrometric analysis. This makes the steel holder not feasible for automated procedures and limits the sensitivity of the analysis in general since only little sample solution can be applied on each spot.

Steel targets have the disadvantage that there is no on target washing for removing contaminating and signal suppressing salts possible.

30

Graphite targets cannot be regenerated, because of the strong absorption of the sample on the surface, which is also true for porous silicone targets. Besides this pure graphite targets have the disadvantage that the mass

- 2 -

spectrometer can be easily contaminated by conductive graphite dust, which will lead to a breakdown of the turbo pumps or the electronic and therefore damage the instrument seriously. Same is true for the so-called liquid matrix, some graphite dispersed in a viscose solvent like glycerol or  
5 silicone oil.

Other holders try to circumvent the problem of sample spreading and small volumes by applying a hydrophobic material on the steel target only leaving a small not coated spot on which the sample concentrates after  
10 evaporation of the solvent (hydrophilic anchors).

A different approach uses a hydrophobic coating and a small spot filled with chromatographic reversed phase C-18 material on a steel target. The samples can then be washed to remove salt after they have been spotted  
15 onto the holder, which increases the quality of the mass spectrometry read-out later. However, because of the nature of the reversed phase material and the strong binding of the peptides to it, it is difficult to regenerate the material after use, which often leads to cross over contamination or loss of the chromatographic material during regeneration  
20 of the target.

Thus, one objective of the invention is to provide a sample holder with a surface to which the peptides bind strongly in order to allow washing of the samples on the holder and which at the same time can be entirely  
25 regenerated leaving no contaminants.

This objective was accomplished according to the invention by providing a sample holder for a mass spectrometer characterized in that it contains a coating comprising silicone and graphite.

30

The inventor has found that applying a thin (e.g. 0.01 to 2mm, ideally 0.2mm) layer comprising a mixture of commercially available silicone with

1 to 70 wt-%, ideally 10- 30 wt-% graphite solves all the above-mentioned problems. The amount of graphite as well as the thickness of the coating can be adjusted according to the respective sample to be measured.

5 The silicone-graphite mix strongly binds peptides, allowing washes, increasing sensitivity, retaining the resolution and feasibility for automation, and is easily removed using conventional silicone removers.

10 Due to the good binding characteristics of the coating according to the invention, the sample can be applied to a smaller surface which leads to a higher concentration of the sample/surface sample holder and thus leads to better results in the mass spectrometry.

15 Any silicone can be used as a silicone component. Preferably, a silicone, which is commercially available, is used. Suitable silicones include any compounds in which Si atoms are connected to O atoms to form chain or net like structures and any remaining valences of Si are connected to hydrocarbon groups. Suitable hydrocarbon groups include C<sub>1</sub>-C<sub>8</sub> alkyl groups, e.g. methyl, ethyl or propyl, C<sub>2</sub>-C<sub>8</sub> alkenyl groups or C<sub>4</sub>-C<sub>15</sub> aryl groups, e.g. phenyl. The hydrocarbon groups preferably contain 1 to 15, in particular 1 to 8 C-atoms and may contain one or more heteroatoms, e.g. selected from N, O or S.

25 The hydrocarbon groups may further comprise substituents, e.g. OH, NH<sub>2</sub>, NO<sub>2</sub>, COOH, C<sub>1</sub>-C<sub>4</sub> alkoxy, halogens or COOR, with R being a C<sub>1</sub>-C<sub>8</sub> hydrocarbon group.

A graphite powder is preferably used as graphite.

30 The manufacture of the coating according to the invention can be effected by mixing a silicone with graphite. This mixture can then be applied to a sample carrier. Preferably, monomers or prepolymers, which can react to a

silicone, are first mixed with graphite, this mixture is applied to a sample carrier and then polymerized on the sample carrier. It is further possible to mix a sample with the graphite and/or silicone components and apply this sample/graphite/silicone mixture as coating to a sample holder. It is possible to add an additional matrix compound to enhance the MS performance, but spectra can also be obtained without the use of such additional matrix substances (cf. Fig. 8).

The sample holder itself can be made of any type of material, preferably of steel. Furthermore, the invention concerns the use of a mixture of silicone and graphite for coating a sample holder for a mass spectrometer.

Furthermore, the invention concerns a method of analyzing a sample in a mass spectrometer comprising the steps

- (a) providing a sample holder containing a coating comprising silicone and graphite,
- (b) applying the sample onto the sample holder and
- (c) performing a mass spectrometry analysis of the sample.

A special advantage of the sample holder coating according to the invention consists in that the sample holder can be washed in order to remove contaminations from the sample, in particular salt contaminations. This is possible, because the sample strongly adheres to the coating that it is not being washed off and on the other hand since contaminations, especially salts, can be removed due to their water solubility. The washing step is preferably carried out with water or aqueous solutions. The sample carrier and/or the method according to the invention is especially suitable in connection with the determination of biomolecules such as proteins, peptides, nucleic acids, steroids, fatty acids, sugars, small molecules ( $M_w < 1000$  Da), especially of proteins and/or peptides.

For mass spectrometric analysis the sample applied to the sample holder is preferably subjected to a laser desorption step.

5 The invention is further explained by the enclosed Figures and the following examples. Figure 1 shows a steel target with 2  $\mu$ l spotted matrix/sample mix. The diameter of a spot is 2,5 mm (resulting in an area of 4,9 mm<sup>2</sup>).

10 Figure 2 shows a silicone/graphite target according to the invention with 2  $\mu$ l spotted matrix/sample mix. The diameter is 1,9 mm (corresponding to an area of 2,8 mm<sup>2</sup>). Thus, the sample is concentrated on a 40% smaller area, which means that there is no search for a good spot on the target necessary. Firing the laser on the silicone-graphite target produces immediately signals, but in the case of the steel target it is most often  
15 necessary to search for good crystallized spots inside the target spot.

Figure 3 shows in an enlargement the fairly wide distributions of the crystals of Matrix/sample on a steel target.

20 Figure 4 shows the homogeneous crystallization of a sample on a silicone/graphite target.

Figure 5 shows a wash step performed on a silicone/graphite target according to the invention. Because of the relatively high hydrophobicity of  
25 the silicone/graphite coating, it is possible to wash such crystallized spots with large amounts of water. Usually, a drop of about 8  $\mu$ l is set on the Matrix/sample spot. The water spot covers just the crystallization area and does not spread further. Because of the large amounts of water one washing step is sufficient since the contaminating salts are very effectively  
30 dissolved in the washing water.

- 6 -

Figure 6 shows the result of a comparison between steel and silicone/graphite targets of samples deriving from a 2D-gel separation and in gel tryptic digestion is shown below. The steel target without washing gave a positive database search identification result of 29% and after washing of the steel target 26%, but with the silicone/graphite coated target, including on-target washing, the positive identification was 79%. See row 1 to 7 (56 samples). Row number 8 was a control containing only extracts from blank gel, no proteins, therefore no positive identification. A list showing the identified proteins is also given.

Figure 7 shows mass spectra obtained using a silicone/graphite coated target according to the invention or a steel target. A comparison of resolution between the steel and silicone/graphite target showed that there are no significant differences (steel 6500 and slightly better silicone/graphite 8100) and therefore the silicone/graphite target is resolution neutral. Comparison of the intensity showed clearly that the silicone/graphite target (here 3350 total ion counts per second) is in average 4 times more sensitive than the steel target (in this example 445 total ion counts per second), which is very important for the analysis of less abundant proteins.

Figure 8 shows a comparison of signal intensities. With the silicone/graphite matrix it is even possible to acquire spectra without using a matrix as shown in Figure 8 where a mixture of 6 peptides was applied to the silicone/graphite coated target without additional matrix.

Figure 9 shows that there are no background signals deriving from the silicone/graphite coating, which is an additional advantage.

ExamplesExample 15     **Materials :**

Silicone: "Knauf kitchen silicone acetat crosslinking (Knauf Bauprodukte GmbH, ID-Nr. 7949)

Graphite: Merck, Graphite, pulver, pure, Order No. 1.04206.2500

10     1.58 g silicone were thoroughly mixed with 0.185 g graphite and immediately transferred onto the MS steel target with the help of a spatula. The target was then pushed through a coating apparatus that produced a defined height of the silicone/graphite layer of 0.2 mm. The polymerisation process took over night.

15

With the help of a sample holder manufactured in that way, the results, which are shown in the Figures, were obtained and compared with a conventional steel target.

20




### Claims

1. A sample holder for a mass spectrometer characterized in that it contains a coating comprising silicone and graphite.
2. Sample holder according to claim 1, wherein the coating layer having a thickness of 0,01 to 1 mm.
3. Sample holder according to claim 1 or 2, wherein the coating contains 1 parts by weight of graphite: 100 parts by weight silicone to 70 parts by weight graphite: 100 parts by weight silicone.
4. Sample holder according to claim 3, wherein the coating contains 10 parts by weight graphite: 100 parts by weight silicone to 30 parts by weight graphite: 100 parts by weight silicone.
5. Sample holder according to any of the preceeding claims, wherein the coating is prepared by polymerizing a silicone-forming monomer or prepolymer in the presence of graphite.
6. Use of silicone and graphite for coating a sample holder for a mass spectrometer.
7. A method of analyzing a sample in a mass spectrometer comprising the steps
  - (a) providing a sample holder containing a coating comprising silicone and graphite,
  - (b) applying the sample onto the sample holder and
  - (c) performing a mass spectrometry analysis of the sample.
8. The method according to claim 7 further comprising a step

(b1) washing the sample holder after application of the sample thereto to remove contaminations from the sample.

- 5        9.     A method according to claim 8, wherein in the washing step b1, salt contaminations are removed.
- 10      10.    The method according to any of claim 7-9, wherein the mass spectrometry analysis comprises a laser desorption step.
- 10      11.    Use of a mixture comprising silicone and graphite in mass spectrometric analysis.
- 15      12.    Use according to claim 6 or 11, wherein the mixture comprises 1 parts by weight of graphite: 100 parts by weight silicone to 70 parts by weight of graphite: 100 parts by weight silicone.
13.    Use according to claim 6 or 11 to 12, wherein the mixture is provided in a thickness of from 0,01 to 1mm.

**Figure 1**

**Figure 2**







Figure 3



**Figure 4**


**Figure 5**






# Figure 6


Steel Target



Steel Target with on-target washing



Silicone/graphite with on-target washing





Sample		Silicone/Graphite	Steel without washing	Steel with washing
1	identified Protein	DNAK_ECOLI (P04475)	data quality too poor	DNAK_ECOLI (P04475)
	score	44,3	-	69,88
	no. of matched peptides	29	-	24
	peptide coverage in %	51	-	45
2	identified Protein	RS1_ECOLI (P02349)	CAD26622	RS1_ECOLI (P02349)
	score	32,25	13,05	21,19
	no. of matched peptides	20	16	19
	peptide coverage in %	36	43	34
3	identified Protein	CH60_ECOLI (P06139)	Q8VTM5	CH60_ECOLI (P06139)
	score	37,25	3,46	-7,18
	no. of matched peptides	27	11	13
	peptide coverage in %	59	25	37
4	identified Protein	TIG_ECOLI (P22257)	Q9HU18	Q9EL49
	score	36,76	-10,13	-8,75
	no. of matched peptides	25	6	8
	peptide coverage in %	62	19	34
5	identified Protein	Q9PN60	Q9F260	Q9ECQ1
	score	2,73	-15,04	3,57
	no. of matched peptides	13	7	14
	peptide coverage in %	34	34	48
6	identified Protein	Q9LYL6	ASNA_ECOLI (P00963)	data quality too poor
	score	3,32	-23,98	-
	no. of matched peptides	15	4	-
	peptide coverage in %	47	18	-
7	identified Protein	IMDH_ECOLI (P06981)	data quality too poor	Q8VCN2
	score	27,05	-	-18,03
	no. of matched peptides	16	-	4
	peptide coverage in %	37	-	12

	Identified Protein	DPPA_ECOLI (P23847)	Q9UID7	DPPA_ECOLI (P23847)
8	score	20,17	-3,91	-7,7
	no. of matched peptides	18	12	18
	peptide coverage in %	33	43	35
	identified Protein	ATPA_ECOLI (P00822)	ATPA_ECOLI (P00822)	ATPA_ECOLI (P00822)
9	score	10,42	35,83	4,47
	no. of matched peptides	10	21	13
	peptide coverage in %	22	44	27
	identified Protein	DLDH_ECOLI (P00391)	Q9EED5	DLDH_ECOLI (P00391)
10	score	25,16	-10,1	42,23
	no. of matched peptides	18	9	15
	peptide coverage in %	46	26	41
	identified Protein	LEU2_ECOLI (P30127)	O64862	AAH27312
11	score	12,95	6,5	-10,88
	no. of matched peptides	14	15	6
	peptide coverage in %	32	36	21
	identified Protein	GLYA_ECOLI (P00477)	GLYA_ECOLI (P00477)	Q9NBV1
12	score	25,62	14,71	-14,19
	no. of matched peptides	13	18	4
	peptide coverage in %	31	41	22
	identified Protein	GLYA_ECOLI (P00477)	Q9EI10	Q78557
13	score	30,02	18,72	-19,77
	no. of matched peptides	19	16	4
	peptide coverage in %	47	51	46
	identified Protein	ENO_ECOLI (P08324)	ENO_ECOLI (P08324)	ENO_ECOLI (P08324)
14	score	44,01	35,81	39,23
	no. of matched peptides	20	20	15
	peptide coverage in %	57	52	30

15	identified Protein	Q9NUJ5	EFTU_ECOLI (P02990)	Q97GN2
	score	25,39	56,52	10,07
	no. of matched peptides	16	14	9
	peptide coverage in %	59	49	12
16	identified Protein	PGK_ECOLI (P11665)	PGK_ECOLI (P11665)	PGK_ECOLI (P11665)
	score	41,74	71,27	0,52
	no. of matched peptides	19	24	17
	peptide coverage in %	61	66	47
17	identified Protein	LIVJ_ECOLI (P02917)	LIVJ_ECOLI (P02917)	LIVJ_ECOLI (P02917)
	score	32,87	36,71	24,79
	no. of matched peptides	11	12	11
	peptide coverage in %	26	27	26
18	identified Protein	EFTS_ECOLI (P02997)	EFTS_ECOLI (P02997)	EFTS_ECOLI (P02997)
	score	21,65	24,42	38,5
	no. of matched peptides	14	13	17
	peptide coverage in %	61	35	50
19	identified Protein	MDH_ECOLI (P06994)	MDH_ECOLI (P06994)	YDCM_ECOLI (P76102)
	score	32,17	14,51	19,45
	no. of matched peptides	16	13	19
	peptide coverage in %	57	45	33
20	identified Protein	CYSK_ECOLI (P11096)	CYSK_ECOLI (P11096)	CYSK_ECOLI (P11096)
	score	39,06	-3,7	17,71
	no. of matched peptides	15	10	13
	peptide coverage in %	57	39	54
21	identified Protein	ZNUA_ECOLI (P39172)	RF1_MYCPU (Q98RA5)	Q9NZE9
	score	32,39	-10,23	28,89
	no. of matched peptides	13	5	13
	peptide coverage in %	48	11	33

22	identified Protein	HISJ_ECOLI (P39182)	HISJ_ECOLI (P39182)	HISJ_ECOLI (P39182)
	score	30,31	-1,62	-12,84
	no. of matched peptides	14	11	11
	peptide coverage in %	50	40	53
23	identified Protein	AROG_ECOLI (P00886)	GLTD_ECOLI (P09832)	AROG_ECOLI (P00886)
	score	21,39	1,42	-1,21
	no. of matched peptides	17	12	10
	peptide coverage in %	50	28	34
24	identified Protein	G3P1_ECOLI (P06977)	Q9VR05	Q9UWZ0
	score	-5,39	-5,89	-20,32
	no. of matched peptides	11	11	4
	peptide coverage in %	42	34	19
25	identified Protein	G3P1_ECOLI (P06977)	G3P1_ECOLI (P06977)	G3P1_ECOLI (P06977)
	score	26,12	-1,71	18,13
	no. of matched peptides	15	10	10
	peptide coverage in %	59	50	39
26	identified Protein	SURA_ECOLI (P21202)	GATM_PIG (P50441)	SURA_ECOLI (P21202)
	score	23,11	31,54	-3,38
	no. of matched peptides	19	12	14
	peptide coverage in %	38	27	32
27	identified Protein	AAH03774	AAL78318	AAH25541
	score	4,16	1,49	-21,36
	no. of matched peptides	15	12	4
	peptide coverage in %	44	47	10
28	identified Protein	TRPB_ECOLI (P00932)	O98097	TRPB_ECOLI (P00932)
	score	14,19	-16,84	-12,15
	no. of matched peptides	14	4	10
	peptide coverage in %	52	43	24

29	identified Protein	TRPB_ECO57 (Q8X7B6)	O31310	Q947A5
	score	10,62	-6,97	-8,87
	no. of matched peptides	16	8	10
	peptide coverage in %	47	35	35
30	identified Protein	CARA_ECOLI (P00907)	Q23796	Q9L357
	score	24,98	-24,67	-7,41
	no. of matched peptides	14	4	8
	peptide coverage in %	43	13	15
31	identified Protein	SUCD_ECOLI (P07459)	Q9EKC6	GALS_ECOLI (P25748)
	score	0,43	-24,49	-2,51
	no. of matched peptides	7	4	9
	peptide coverage in %	30	22	31
32	identified Protein	SUCD_ECOLI (P07459)	Q92WT8	SUCD_ECOLI (P07459)
	score	11,13	-16,84	-3,03
	no. of matched peptides	11	5	10
	peptide coverage in %	44	18	45
33	identified Protein	Q14533 (Keratin)	K1M2_SHEEP (P25690)	Q28582 (Keratin)
	score	19,46	11,95	38,01
	no. of matched peptides	21	14	20
	peptide coverage in %	40	31	38
34	identified Protein	RS2_ECOLI (P02351)	Q9XW07	RS2_ECOLI (P02351)
	score	20,9	-4,04	10,54
	no. of matched peptides	12	11	8
	peptide coverage in %	38	40	38
35	identified Protein	ARGT_ECOLI (P09551)	ARGT_ECOLI (P09551)	Q8VXK3
	score	22,18	27,81	8,63
	no. of matched peptides	14	16	9
	peptide coverage in %	50	63	23


	identified Protein	FLIY_ECOLI (P39174)	FLIY_ECOLI (P39174)	DAPB_ECO57 (P58209)
36	score	53,53	21,99	-4,59
	no. of matched peptides	14	9	8
	peptide coverage in %	49	27	47
	identified Protein	Q9RCJ3	Q9K2N9	Q9ZTP0
37	score	20,14	-7,26	0,31
	no. of matched peptides	15	9	11
	peptide coverage in %	35	22	23
	identified Protein	TPIS_ECOLI (P04790)	Q33364	O44593
38	score	19,47	3,07	6,67
	no. of matched peptides	12	10	14
	peptide coverage in %	41	23	26
	identified Protein	Q9TQ74	Q9EF71	ARTJ_ECOLI (P30860)
39	score	10,09	-19,47	-9,26
	no. of matched peptides	11	4	8
	peptide coverage in %	28	14	45
	identified Protein	TRPA_ECOLI (P00928)	Q9LV79	TRPA_ECOLI (P00928)
40	score	7,35	-12,23	9,38
	no. of matched peptides	10	6	13
	peptide coverage in %	46	34	56
	identified Protein	AAL90243	Q67371	Q20218
41	score	1,72	-2,01	16,52
	no. of matched peptides	12	7	16
	peptide coverage in %	37	29	47
	identified Protein	PMG1_ECOLI (P31217)	PMG1_ECOLI (P31217)	PMG1_ECOLI (P31217)
42	score	27,99	0,5	6,54
	no. of matched peptides	15	7	10
	peptide coverage in %	56	27	37

43	identified Protein	YODA_ECOLI (P76344)	Q98NV8	YODA_ECOLI (P76344)
	score	7,16	4,63	3,93
	no. of matched peptides	5	10	5
	peptide coverage in %	21	25	19
44	identified Protein	Q9E6M0	YNEC_ECOLI (P76144)	Q9F753
	score	26,85	-22,18	-3,97
	no. of matched peptides	13	4	7
	peptide coverage in %	34	13	18
45	identified Protein	GLNH_ECOLI (P10344)	GLNH_ECOLI (P10344)	Q90EU5
	score	17,29	-3,65	-17,39
	no. of matched peptides	10	8	4
	peptide coverage in %	50	32	35
46	identified Protein	AHPC_ECOLI (P26427)	AHPC_ECOLI (P26427)	Q8XKG2
	score	42,25	10,75	-19,82
	no. of matched peptides	14	10	4
	peptide coverage in %	73	51	15
47	identified Protein	Q9NV51	BAB88645	O92108
	score	-5,01	-13,79	-19,45
	no. of matched peptides	4	8	4
	peptide coverage in %	8	26	41
48	identified Protein	Q96M01	DNAJ_RHOSP (O08356)	Q9X9N4
	score	-1,18	-1,91	-13,11
	no. of matched peptides	9	9	4
	peptide coverage in %	23	22	11
49	identified Protein	O01622	Q8U7H2	O37399
	score	0,46	-4,45	-12,03
	no. of matched peptides	7	10	5
	peptide coverage in %	16	40	9

	identified Protein	Q90ZT5	Q74651	Q9IZ27
50	score	-5,46	-9,44	-11,88
	no. of matched peptides	9	11	6
	peptide coverage in %	34	38	21
	identified Protein	RL9_ECOLI (P02418)	Q97M24	RL9_ECOLI (P02418)
51	score	2,95	-2,2	2,1
	no. of matched peptides	6	8	7
	peptide coverage in %	38	34	37
	identified Protein	METE_ECOLI (P25665)	GSHB_BUCAI (P57612)	YDIP_ECOLI (P77402)
52	score	14,54	-2,18	-20,68
	no. of matched peptides	26	10	4
	peptide coverage in %	33	33	10
	identified Protein	METE_ECOLI (P25665)	METE_ECOLI (P25665)	O95958
53	score	25	7,59	8,77
	no. of matched peptides	25	14	18
	peptide coverage in %	34	15	40
	identified Protein	METE_ECOLI (P25665)	O45505	Q9HVV0
54	score	44,23	-18,87	-14,26
	no. of matched peptides	26	4	5
	peptide coverage in %	40	12	17
	identified Protein	EFG_ECOLI (P02996)	EFG_ECOLI (P02996)	Q9EHF6
55	score	42,63	4,42	-4,51
	no. of matched peptides	25	19	9
	peptide coverage in %	45	33	13
	identified Protein	Q9X788	data quality too poor	Z265_MOUSE (Q9R020)
56	score	3	-	-0,92
	no. of matched peptides	10	-	7
	peptide coverage in %	41	0,92	12
Number of positive identified proteins		43	19	25
Percentage of positive identified proteins		77	34	45
Total Time for data acquisition/database search		145.6 min	257.6 min	257.6 min
acquisition/database search time per sample		2.6 min	4.6 min	4.6 min




Figure 7



# Figure 8

Comparison of the signal intensities

silicone/graphite coated target



steel target




Figure 9

